

Areterna LLC



# THE STORY OF

# mRNA


<b>GMP</b>	<b>CHEMICAL</b>	<b>PLASMID</b>	<b>SEQUENCE</b>
	<b>IVT</b>	<b>LNP</b>	<b>QUALITY CMC</b>




# Background

The success of Pfizer and Moderna COVID vaccines has catapulted mRNA to the forefront of drug development. The discovery of messenger RNA (mRNA) dated back to 1961, however, the molecule is not a viable therapeutic modality until the discovery of modified uridine and lipid nanoparticle (LNP) encapsulation. Both modified uridine and Cap1 structure help mRNA evade cell's innate immune system and thereby improve mRNA's stability. LNP protects mRNA in circulation and aids its cellular uptake and endosomal escape. Compared to other biologics, mRNA offers the benefits of simple manufacturing, quick design and no risk of genome integration. Unlike AAV or antibody drugs, the manufacturing of mRNA does not involve any cells. It starts with in-vitro transcription of a DNA template, followed by purification, LNP encapsulation and sterile fill-finish. As a new therapeutic modality, mRNA offers great potential to treat and prevent diseases, from vaccine, to cancer treatment, to protein replacement, gene editing, cell therapy and many more.

## Manufacturing Process of mRNA Vaccine



# What Areterna Can Do



## Products

- Cap Analogs
- NTPs and Modified NTPs
- GMP and Non-GMP Grade
- Custom Options

- mRNA Synthesis Kits
- Capped/Uncapped mRNA Standards
- Residual dsRNA Detection Kits
- Impurity Standards
- Stocked Tool mRNAs
- LNP Screening Plate





## Services

- mRNA Sequence Design
- mRNA Custom Synthesis
- IVT Process Development
- mRNA Testing
- Analytical Development
- LNP Formulation Screening
- LNP Encapsulation
- IND Filing Support

## Technology Platforms

- Sequence Design Platform
- Nucleotide Chemistry Discovery Platform
- IVT Manufacturing Platform
- LNP Encapsulation Platform





## 1 Sequence Design

**Services:** Target gene screening  
Sequence design and optimization

## 2 DNA Template Construction

**Products:** dNTPS  
**Services:** Plasmid construction and linearization

## 3 In-vitro Transcription


**Products:** Co-transcription kits  
**Services:** Custom mRNA synthesis

## 5 Analytical Testing

**Products:**  
Impurity standards  
Cap1&Uncap standards  
Residue test kits  
**Services:** mRNA release test  
mRNA stability testing  
Assay development and trans (e.g. cleavage probe design)

## 4 LNP Encapsulation

**Products:** Stocked Tool mRNAs  
**Services:** LNP formulation screening  
LNP encapsulation service



## IND Support

9

**Services:** Preparation of regulatory documents  
Consultation on process compliance

7

## Process Development

**Services:** IVT Process development  
Analytical method development and validation

8

## GMP Manufacturing

**Products:** Cap analogs  
Cap analog intermediate  
Modified nucleotides  
Nucleotides  
Enzymes

6

## Function Testing

**Products:** Cell transfection kit



### **Areterna Services**

Target Gene Screening



### **Areterna Services**

Sequence Design And Optimization

# **Chapter One**


## **Sequence Design**

mRNA is more than a coding sequence. From 5' end, there is a cap structure that binds eIF4E to initiate protein translation, then 5' UTR for ribosome scanning, followed by coding sequence, then 3' UTR that regulates protein expression, and poly-A tail to stabilize mRNA. Each segment can be optimized for optimal protein expression, reduced immunogenicity and better manufacturability. Selecting the right gene or protein to modulate is also critical. One of the challenges in cancer immunotherapy is to identify the appropriate antigen that elicits strong and durable immune response.

# Target Gene Identification


# ARETERNA SERVICES

Target gene identification is the starting point for an mRNA drug or vaccine and can be a complex process. Poor choice of target gene can lead to failed drugs and wasted resources. The diagram below is an example of tumor neoantigen selection process. Areterna can help customers to screen and identify the right targets.



# Sequence Design and Optimization

Areterna has its own proprietary algorithm for sequence design and optimization. All we need is the amino acid sequence, we can add UTRs, cap and poly-A, and adapt the sequence to specific species if veterinary vaccine is the intended use. We use our proprietary cap analogs to introduce Cap1 structure. The approaches to optimizing UTRs and CDS are listed below-



**Areterna Products**  
dNTPs



**Areterna Services**  
Plasmid Construction And Linearization



## Chapter Two


### DNA Template Construction

To synthesize mRNA, a DNA template is needed. It can be a linearized plasmid or a PCR product. For small scale synthesis, PCR products are often used. However, PCR is not scalable, so plasmid becomes the template of choice for large scale GMP manufacturing. Regardless of DNA source, the template should contain T7 Polymerase promoter sequence in addition to the 5 parts of mRNA. One confusing point for mRNA novices is that the cap is not included in the DNA template as the inverted G is added by the cap analog or enzymatic reaction. If linearized plasmid is used, care should be taken to avoid linearization site in the sequence.



Areterna offers high quality dNTPs.

- ✓ High purity | Purity verified by HPLC, higher than 99%
- ✓ High quality | Ensuring batch stability and freeze-thaw stability
- ✓ High sensitivity | No DNase, RNase and cleavage enzyme residues
- ✓ High specificity | No bacterial or human genome contamination





## dNTPs Product List

Cat. No.	Product Name	Cat. No.	Product Name
dNTP001	HiPure® dATP	dNTP009	HiPure®Cy3-dUTP
dNTP002	HiPure® dGTP	dNTP010	HiPure®Biotin-11-dUTP
dNTP003	HiPure® dCTP	dNTP011	HiPure® 3'-ONH <sub>2</sub> -dATP
dNTP004	HiPure® dTTP	dNTP012	HiPure® 3'-ONH <sub>2</sub> -dGTP
dNTP005	HiPure® dUTP	dNTP013	HiPure® 3'-ONH <sub>2</sub> -dCTP
dNTP006	HiPure® dNTP Mix (25mM each)	dNTP014	HiPure® 3'-ONH <sub>2</sub> -dTTP
dNTP007	HiPure® dNTP Mix (10mM each)	dNTP015	HiPure® Fluorescein-12-dUTP
dNTP008	HiPure® dNTP Mix (2.5mM each)		

## Plasmid Construction and Linearization

To help customers with their mRNA drug development, Areterna offers plasmid synthesis service from microgram (ug) to milligram (mg) quantities. The linearized plasmid will be sequence verified and ready to go into IVT reaction.

### Flow Chart of Plasmid Construction and Linearization



# Chapter Three

## In Vitro Transcription

In-Vitro Transcription (IVT) is the core of mRNA synthesis. It utilizes T7 or SP6 RNA Polymerase to incorporate NTPs into the growing mRNA chain. Cap can be added enzymatically post-IVT or in a one-pot co-capping reaction. The enzymatic capping is a 4-step process using vaccinia capping enzyme and 2'-O-methyl transferase, the capping efficiency can be very high if the process is optimized. One-pot co-capping requires less steps and less unit of operation. The capping efficiency varies depending on the cap analogs used. Anti-Reverse-Capping-Analog (ARCA) is the first-generation cap analog which yields a Cap0 structure with moderate capping efficiency. GpppAG cap analog and its variation such as Areterna's GAG (UNA) and GAG (ENE) offer high capping efficiency and Cap1 structure, therefore have become the cap analogs of choice for many applications. To reduce the immunogenicity of mRNA, modified uridine or cytidine can be used in the IVT reaction.


**Areterna Products**  
Mrna Co-transcription Kits



**Areterna Services**  
mRNA Custom Synthesis



In order to help bench scientists to quickly synthesize mRNA for research, Areterna introduced a T7 Co-transcription mRNA Synthesis Kit, with our novel Cap1 cap analogs and N1-methyl pseudo uridine. The kit contains all components needed for IVT except the DNA template.




- ✓ **Ease of use** | Cap1 mRNA can be synthesized in one step
- ✓ **High yield** | up to 200 µg of mRNA can be generated from 1µg DNA template
- ✓ **High capping efficiency** | capping efficiency of > 95% can be achieved
- ✓ **Reduced immunogenicity** | by replacing wild-type UTP with N1-Me-pUTP

### Recommended IVT Conditions


10*Reaction Buffer	2 µL
100mM ATP Solution	2 µL
100mM CTP Solution	2 µL
100mM GTP Solution	2 µL
Reagent A	2 µL
DNA Template	1 µg
Reagent B	2 µL
T7 RNA polymerase (250U/µL)	1 µL
Inorganic Pyrophosphatase (1U/µL)	0.04 µL
RNase inhibitor (40U/µL)	1 µL
RNase free ddH <sub>2</sub> O	up to 20 µL

## Co-transcription Kit Performance


The HiSynth® T7 Co-transcription mRNA Synthesis Kit works with short and long DNA templates. We've seen high protein expression from mRNA constructs ranging from 1.4kk to 10kb. Compared to enzymatic capping method, the kit yields the similar amount of mRNA which shows higher activity in the data below.




Accommodate Templates with various sizes



Similar Yiled as Enzymatic Capping



High activity of mRNA products



Capping efficiency  $\geq 95\%$


## The IVT and Capping Kit Product List

Cat. No.	Product Name
10110U	HiSynth® T7 High Yield RNA Synthesis Kit (AGCU)
10110N	HiSynth® T7 High Yield RNA Synthesis Kit (AGCN)
10111-4011	HiSynth® T7 Co-transcription RNA Synthesis Kit(4011)
10111-5011	HiSynth® T7 Co-transcription RNA Synthesis Kit(5011)
10111-6011	HiSynth® T7 Co-transcription RNA Synthesis Kit(6011)
10111-7011	HiSynth® T7 Co-transcription RNA Synthesis Kit(7011)
10111-8011	HiSynth® T7 Co-transcription RNA Synthesis Kit(8011)
10111-9011	HiSynth® T7 Co-transcription RNA Synthesis Kit(9011)
10112	HiSynth® Cap 1 Capping System

# 3 ARETERNA SERVICES


## mRNA Custom Synthesis

For customers who do not have time or expertise to synthesize mRNA, Areterna offers custom synthesis service using our HiTrans® IVT platform. We can start with digital sequence and deliver purified mRNA in milligram (mg) and gram (g) quantities with accompanying QC data.




### Examples of Custom Synthesized mRNA

Quality of Linearized Plasmid




Capping Efficiency by LC-MS



mRNA Integrity by Fragment Analyzer




Poly A tail length and distribution by LC-MS



# Chapter Four

## LNP Encapsulation

Naked mRNA is vulnerable to digestion by circulating nucleases and has poor intracellular uptake properties due to its size and charge. Lipid nanoparticle (LNP) is used to protect mRNA and facilitate its delivery. LNPs are typically composed of four lipids - cationic or ionizable lipids, cholesterol, phospholipids, and PEG-lipids. Ionizable lipids are essential for mRNA complexation by forming a core around negatively charged mRNA. They remain neutral in circulation but become protonated in endosome which facilitate mRNA's endosomal escape. PEG-lipids play an important role in extending LNP's half-life in circulation and preventing particle aggregation. Phospholipids and cholesterol contribute to the structural integrity of the LNP. The generation of LNP relies on rapid mixing of mRNA in aqueous phase with lipids dissolved in ethanol.



**Areterna Products**  
Stocked Tool mRNAs



**Areterna Services**  
LNP Formulation Screening



**Areterna Services**  
LNP Encapsulation




Areterna stocks a list of “Tool” mRNAs for LNP encapsulation and delivery research, including eGFP mRNA, FLuc mRNA, OVA mRNA and Fluc-eGFP mRNA. The beauty of Fluc-eGFP mRNA is that it can be assayed both in cells (eGFP) and in mice (FLuc). The construct is longer than FLuc and can be a better surrogate for longer mRNAs.



### FLuc-eGFP mRNA Data

01



#### mRNA Integrity



Fluc-eGFP mRNA integrity > 90%  
by Fragment Analyzer

02

#### Capping Efficiency




Capping Efficiency > 95% by LC-MS



03

## mRNA Activity



Cells were transfected with Fluc-eGFP mRNA. Both GFP and luciferase show strong signal 24 hours post transfection.

## Stocked Tool mRNAs

### PRODUCT NAME

HiSignal® eGFP mRNA

HiSignal® Fluc mRNA

HiSignal® Fluc-eGFP mRNA

HiSignal® mCherry mRNA mRNA


HiSignal® OVA mRNA

# LNP Formulation Screening LNP Encapsulation Service

Areterna has a library of ionizable lipids and makes it available to customers for screening. The LNP formulation screening can be conducted at Areterna or at customer's site. To help customers with mRNA research, Areterna also offers LNP encapsulation service should customers choose to do so.

## Schematic Diagram of Typical LNP Formula

LNPs are typically composed of four lipids - cationic or ionizable lipids, cholesterol, phospholipids, and PEG-lipids.




Reference: Acta Biomater. 2021 Sep 1;131:16-40. doi: 10.1016/j.actbio.2021.06.023.

## LNP Encapsulation


Particle Size 60-150nm

PDI < 0.2

Encapsulation Efficiency > 90%



Particle Size and PDI



Zeta Potential

# Chapter Five

## Analytical Testing


To release an mRNA drug substance or drug product, a set of testing needs to be performed to monitor critical quality attributes, process-related residuals and safety. Some tests are standard, such as appearance, pH, endotoxin, bioburden and sterility. Key attributes that are unique to mRNA include characterization of the 5' cap and poly-A tail to ensure mRNA translation in vivo, as well as assessment of dsRNA as a potential cause for immunogenicity. If the mRNA is encapsulated with LNP, additional tests need to be done to measure particle size, polydispersity index, %encapsulation, lipid identify, purity and content.




**Areterna Products**  
Impurity Standards




**Areterna Services**  
mRNA Analytical Testing




**Areterna Products**  
Cap1&uncap Standards



**Areterna Services**  
LNP Analytical Testing




**Areterna Products**  
Residue Test Kits



As customers progress from early clinical phase to commercial phase, quality control of starting materials becomes important. Critical impurities need to be identified and tested. To help customers validate impurity testing methods, Areterna offers a variety of impurity standards.


### Cap Analog Impurity Standards



### Paired Capped /Uncapped mRNA Standards

To assay capping efficiency, a cleavage probe is needed to bind to the 5'UTR of the mRNA and cleave off a short sequence so that the presence or absence of one guanosine can be resolved by LC-MS. Areterna can design the cleavage probe for customers, or develop the capping assay and tech transfer the method, or even run the capping assay for customers. We also custom synthesize paired capped & uncapped mRNA standards for customers to validate their capping assays.

### Paired Capped /Uncapped mRNA Standards



# ELISA Kits for Detecting IVT Process Residues

Process-related impurities need to be monitored, particularly dsRNA which is considered as a CQA as it can trigger innate immune response. Other IVT related impurities include T7 RNA Polymerase, DNase I, RNAase Inhibitor, Inorganic Pyrophosphatase. Areterna offers ELISA kits to test these residues.

## ELISA Kits for Detecting IVT Process Residues



Product Name	Intended Use
T7 RNA Polymerase Test Kit	Quantitative detection of T7 RNA Polymerase content in mRNA stock solution.
DNase I Test Kit	Quantitative detection of DNase I content in mRNA stock solution.
RNase Inhibitor Test Kit	Quantitative detection of RNase Inhibitor content in mRNA stock solution.
Pyrophosphatase Inorganic Test Kit	Quantitative detection of Pyrophosphatase Inorganic content in mRNA stock solution.
dsRNA Test Kit	Quantitative detection of dsRNA content in mRNA stock solution.

### mRNA Vaccines/Drugs Enzyme Residue Test Product List

Cat. No.	Product Name
T70001	T7 RNA Polymerase Test Kit
DI0001	DNase I Test Kit
RI0001	RNase Inhibitor Test Kit
PI0001	Pyrophosphatase Inorganic Test Kit
DS0001	dsRNA Test Kit

# 5 ARETERNA SERVICES

## Analytical Method Development mRNA Release Testing and Stability Testing

Areterna provides full panel release testing and stability testing for mRNA and mRNA/LNP. We can also develop analytical methods such as capping and tailing assays, and tech transfer to customers.


### mRNA QC Panel

Test	Method
Appearance	Visual inspection
pH	USP <791>
Bioburden	USP <61>
Endotoxin	USP <85>
Concentration	UV Spectrophotometry (A260)
Identity	Sequencing (Outsource)
Purity & Integrity	CGE /HPLC
Capping Efficiency	LC-MS
PolyA Tail Length	LC-MS
Residual pDNA	qPCR
Residual Protein	NanoOrange /ELISA
Residual NTPs and Cap Analog	HPLC
dsRNA	ELISA

### mRNA/LNP QC Panel

Test	Method
Appearance	Visual inspection
Osmolality	USP <785>
pH	USP <791>
Particulate Matter	USP <788>
Endotoxin	USP <85>
Lipid Identity & Content	HPLC- CAD
RNA Content and % Encapsulation)	RiboGreen Assay
Polydispersity Index & Particle Size	Zetasizer
Residual Solvents	GC (FID)
Biological Activity	Flow /ELISA
Sterility	USP <71>

### mRNA/LNP Analytical Platforms



LC-MS



FA5200



HPLC



LC-MS



Tecan Spark



Zetasizer Lab



**Areterna Products**  
Cell Transfection Kit



## Chapter Six

### Function Testing

During the process of mRNA drug product development, cell based assays are often performed to quickly verify the protein expression from mRNA. Different mRNA constructs will be tested in cells for their activity, potency and toxicity before advancing the final candidates to animal testing.


HiTrans® mRNA Easy cell transfection kit relies on LNP to complex mRNA and deliver it into cells. The kit has been tested on many commonly used cell lines and was found to offer-

- ✓ High transfection efficiency
- ✓ Low toxicity
- ✓ Good repeatability
- ✓ Cost-effectiveness




### Performance Data

Transfection of eGFP mRNA  
using HiTrans® mRNA Easy kit



Transfection of eGFP mRNA  
using commercial transfection kit



HeLa cells on a 48-well plate were transfected with the HiTrans® mRNA Easy kit or a commercial transfection kit. The eGFP expression was analyzed 24h after transfection. Cells transfected with HiTrans® mRNA Easy kit much stronger eGFP expression.





## Areterna Services

IVT Process Development



## Areterna Services

Analytical Method Development And Validation




# Chapter Seven

## Process Development

To advance the research-lab-made mRNA to the clinical stage, one needs to consider the manufacturability of the construct, batch-to-batch consistency and scalability of the process. Critical Quality Attributes (CQAs) and Critical Process Parameters (CPPs) need to be identified, and the non-compendial analytical methods need to be qualified and validated. Sometimes the construct may need to be redesigned to fit the cGMP manufacturing. Regardless, the IVT process and downstream purification process need to be investigated and optimized to deliver consistent CQAs.

There is no IVT process that fits all. Sometimes trade-offs need to be made to achieve the best purity, highest yield or the lowest dsRNA. Depending on the CQAs desired, many process parameters can be modified, such as [Mg<sup>2+</sup>], [NTP], reaction temperature, amount of DNA template, and amount of T7 polymerase. To scale up from 20ul reaction in microcentrifuge tubes to 10L reaction in bioreactors, one also needs to consider mixing, heating, and reagent addition, etc. Areterna can help customers with process development and scale up.



## Analytical Method Development and Validation

There are a few assays that are unique to mRNA, such as capping efficiency, poly-A tail length and dsRNA content. Capping efficiency assay is construct specific and requires designing of a cleavage probe specific to customer's mRNA. Areterna can develop the analytical methods and transfer to customers.

Test	Method
Residue DNA	qPCR
Residue Protein	Nano Orange
Capping Efficiency	HPLC-MS
mRNA purity	HPLC/CGE
dsRNA content	ELISA
polyA Tail Length	LC-MS
Identity by Oligo Mapping	LC-MS

# Chapter Eight

## GMP Manufacturing

For clinical trials and commercial use, the mRNA needs to be manufactured following cGMP practice to ensure safety and effectiveness of the product. As mRNA is considered as an Advanced Therapy Medicinal Product (ATMP), there is a phase appropriate approach reflecting the increasing degree of control as the drug moves from Phase I to commercial launch. cGMP manufacturing means tighter control of raw materials and requirement for the raw material supplier to have a robust quality system to support high grade products. Areterna offers GMP grade cap analogs, NTPs and modified NTPs, with Drug Master Files on a selected number of products.

**Areterna Products**  
Cap Analogs



**Areterna Products**  
Modified Nucleotides



**Areterna Products**  
Nucleotides




**Areterna Products**  
Enzymes




Cap analog is considered as a critical starting material for the synthesis of mRNA. In cytosol, the cap structure binds to eIF4E to initiate translation. It also helps to protect mRNA from exonuclease and mediate the turnover of mRNA.


Using AI-assisted structure design and molecular docking simulation, Areterna has developed several novel cap analogs that offer comparable capping efficiency and yield as GpppAG.




## Areterna's Cap Analogs




CAP 4




CAP 5




CAP 6



CAP 7




CAP 8




CAP 9

Of the new cap analogs, CAP 4 and CAP 5 have shown prolonged protein expression in cells and mice.



eGFP mRNA



luciferase mRNA


### Cap Analogs Product List

Cat. No.	Product Name	Cat. No.	Product Name
CAP3130	ARCA	CAP7011	CAP 7
CAP4011	CAP 4	CAP8011	CAP 8
CAP5011	CAP 5	CAP9011	CAP 9
CAP6011	CAP 6		

The commercial success of mRNA vaccine is largely attributable to the discovery of modified uridine by Dr. Kariko and Dr. Weissman. The uridine-rich RNA is recognized by toll-like-receptors as virus and triggers immune response. Replacing wild type uridine with pseudo-uridine or N1-methyl pseudo uridine will improve the stability and translational capacity of mRNA. Areterna offers GMP grade pseudo-UTP and N1-Me pUTP\*.


## Product Recommendation

Pseudo-UTP



Cat. No.: PUTP001

Purity 99.54%



- ✓ High purity | >99% by HPLC
- ✓ Low immunogenicity | Low immunogenicity
- ✓ High quality | GMP grade/non-GMP grade
- ✓ High manufacturing capacity | Available in kilogram quantities

## Modified Nucleotides Product List


Cat. No.	Product Name	Cat. No.	Product Name
PUTP001	HiPure® Pseudo-UTP 100mM Sodium solution	F12UTP001	HiPure® Fluorescein-12-UTP
PUTP001T	HiPure® Pseudo-UTP 100mM Tris solution	CY3UTP001	HiPure® Cy3-UTP
PUTP002T	HiPure® Pseudo-UTP 200mM Tris solution	2NHATP001	HiPure® 2'-NH2-ATP
*NMPUTP001	HiPure® N1-Me-pUTP 100mM Sodium solution	MOEATP001	HiPure® 2'-MOE-ATP
*NMPUTP001T	HiPure® N1-Me-pUTP 100mM Tris solution	3M5OMUTP001	HiPure® 3-Me-5-OMe-UTP
*NMPUTP002T	HiPure® N1-Me-pUTP 200mM Tris solution	DM6ATP001	HiPure® 6-N,N-Dimethyl-ATP
N7MGTP001	HiPure® N7-Me-GTP	SATP001	HiPure® ATP(αS)
5MCTP001	HiPure® 5-Me-CTP	DSATP001	HiPure® SPS ATP
M6ATP001	HiPure® m6A-ATP	5AUTP001	HiPure® 5-Aminopropargyl UTP
5OMUTP001	HiPure® 5-OMe-UTP	3OAAATP001	HiPure® 3'O-Allyl-ATP
CY5UTP001	HiPure® Cy5-UTP	3OAGTP001	HiPure® 3'O-Allyl-GTP
		3OACTP001	HiPure® 3'O-Allyl-CTP
		3OAUTP001	HiPure® 3'O-Allyl-UTP

\*N1-Me-pUTP has DMF filed at FDA


Areterna offers GMP grade NTPs in sodium salt or Tris solution, at 100mM or 200mM.

### High Purity Nucleotides ≥99%


ATP | Purity 99.89%




GTP | Purity 99.62%



CTP | Purity 99.53%



UTP | Purity 99.60%



High purity

>99% by HPLC



High quality

GMP grade/ non-GMP grade



High manufacturing capacity

Available in kilogram quantities

### Nucleotides Product List

Cat. No.	Product Name	Cat. No.	Product Name
ATP001	HiPure® ATP 100mM Sodium solution	CTP001T	HiPure® CTP 100mM Tris solution
GTP001	HiPure® GTP 100mM Sodium solution	UTP001T	HiPure® UTP 100mM Tris solution
CTP001	HiPure® CTP 100mM Sodium solution	ATP002T	HiPure® ATP 200mM Tris solution
UTP001	HiPure® UTP 100mM Sodium solution	GTP002T	HiPure® GTP 200mM Tris solution
ATP001T	HiPure® ATP 100mM Tris solution	CTP002T	HiPure® CTP 200mM Tris solution
GTP001T	HiPure® GTP 100mM Tris solution	UTP002T	HiPure® UTP 200mM Tris solution

\*ATP has DMF filed at FDA

\* CTP has DMF filed at FDA

\* GTP has DMF filed at FDA



To provide one-stop shop for customers, Areterna offers enzymes required for mRNA synthesis, from linearization enzymes, Dnase I, to T7 RNA Polymerase and vaccinia capping enzyme.

## Enzymes Product List

Cat. No.	Product Name	Cat. No.	Product Name
10201-0.5K 10201-2.5K 10201-10K	HiTrans® BspQI	1030M-0.1 1030M-1 1030M-10	HiTrans® T7 RNA polymerase mix
10202	HiTrans® Bsal	10221-001 10221-010 10221-025 10221-500	HiPure® Proteinase K
10203	HiTrans® XbaI	10511-005 10511-040 10511-450	HiPure® RNA Isolation beads
10204	HiTrans® XhoI	10501-005 10501-060 10501-450	HiPure® DNA Isolation beads
10301-10K 10301-100K 10301-250K	HiTrans® T7 RNA polymerase	10401-2K 10401-10K 10401-100K 10401-5M	HiTrans® Vaccinia capping Enzyme
10302A-0.5 10302A-1 10302A-5	10×Transcription Buffer	10402-001 10402-010 10402-025 10402-500	10×Capping buffer
10302B-0.5 10302B-1 10302B-5	10×Transcription Buffer	10403-10K 10403-50K 10403-250K 10403-20M	HiTrans® 2-O'-Methyltransferase
10303-10 10303-100 10303-1K 10303-40K	HiTrans® Inorganic Pyrophosphatase		
10304-10K 10304-20K 10304-40K 10304-100K 10304-1M	HiTrans® Murine RNase Inhibitor		
10211-0.5K 10211-2K 10211-10K	HiTrans® DNase I		



### **Areterna Services**

Preparation Of Regulatory Documents



### **Areterna Services**

Consultation On Process Compliance


## **Chapter Nine**

### **IND Support**

IND submissions must be filed with regulatory agencies before human clinical trials can start. The submission includes detailed information on the drug's chemical composition, proposed clinical trial design and safety and efficacy data from preclinical studies. One challenge in IND submissions is the requirement for extensive documentation, such as Chemistry, manufacturing and control (CMC) data.

# Preparation of Regulatory Documents ARETERNA SERVICES

Areterna has worked with many customers in China and assisted them in their IND filing. We can provide data on our raw materials and help to write certain sections for Module 2.3 and Module 3.



PICTURE SOURCE: WWW.ICH.ORG/PAGE/CTD

## Consultation on Process Compliance

### Regulations and Guidelines for mRNA Vaccines and Drugs



World Health Organization  
 Evaluation of the quality, safety and efficacy of messenger RNA vaccines for the prevention of infectious diseases: regulatory considerations

In October 2021, WHO issued guidance on "Evaluation of the Quality, Safety and Efficacy of messenger RNA Vaccines for the Prevention of Infectious Diseases: Regulatory Considerations" which outlined the regulatory considerations for the quality, safety and efficacy of mRNA vaccines.




In April 2023, the US Pharmacopoeia (USP) published the second edition of "Analytical Procedures for mRNA Vaccine Quality". The draft guideline proposed testing panels for plasmid DNA, mRNA drug substance and mRNA drug product with both release assays and characterization assays.



In August 2020, CDE released "Technical Guidelines for Pharmaceutical Research of mRNA Vaccines for Prevention of Novel Coronavirus (Trial Implementation)" where CMC requirements for mRNA drug substance and drug product are outlined, covering manufacturing process, product characterization, specification and stability testing among others.



In May 2022, CDE published "Technical Guidelines for Pharmaceutical Research of ex-Vivo Gene Editing (Trial Implementation)". The guideline covers mRNA as one of the non-viral gene editing methods and provides guidance on mRNA manufacturing, quality attributes and testing for activity and safety.



## Materials and Services to Enable mRNA Vaccines and Therapeutics

**Areterna LLC (A Synthgene Company)**

15810 Gaither Drive, Suite 230, Gaithersburg, MD 20877

[www.areterna.com](http://www.areterna.com) [info4@areterna.com](mailto:info4@areterna.com)



ARETERNA CORPORATION · 15810 Gaither Drive, Suite 230, Gaithersburg, MD 20877

[www.areterna.com](http://www.areterna.com) · 2023 ARETERNA CORPORATION · ALL RIGHTS RESERVED · PRICES AND SPECIFICATIONS SUBJECT TO CHANGE WITHOUT PRIOR NOTICE · PRINTED IN CHINA

Distributed in Europe by

**BIORON GmbH**

In den Rauhweiden 20  
67354 Römerberg  
Germany

Phone +49 6232 298 45-0

Fax +49 6232 29845-29

[info@bion.net](mailto:info@bion.net)

[www.bion.net](http://www.bion.net)

