



Just fine Molecular Biology

Phage Lambda DNA

Description: Phage Lambda d857 *Sam7* DNA is isolated from an infected *E. coli* strain W3350. Restriction enzyme-digested lambda DNA (48,502bp) generates molecular weight size markers routinely used in gel analysis of nucleic acids. The complete nucleotide sequence has been determined (1).

Concentration: 0,2-0,5 mg/ml

Storage Buffer: 10mM Tris-HCl (pH 8,0), 10mM NaCl, 1mM EDTA.

Storage: Store at -20°C.

1. Sanger, F. *et al.* (1982) *J. Mol. Biol.* **162**, 729.

Catalog #	Pack size
301005	0,05 mg
301025	0,25 mg
301100	1 mg

VersionER17.04.08