



Just fine Molecular Biology

Phage Lambda DNA dam-, dcm-

Description Phage Lambda d857 *Sam7* DNA is isolated from an infected *E. coli* strain deficient in adenine and cytosine methylases (dam-dcm-). Restriction endonucleases that are sensitive to dam and/or dcm methylation give complete digests with this DNA.

Storage Buffer 10mM Tris-HCl (pH 7.8), 10mM NaCl, 1mM EDTA.

Storage Store at -20°C.

Loading The recommended amount of size marker to load on an agarose gel is 0,4-0,6 µg/lane

<u>Cat#</u>	<u>Pack size</u>	<u>Price, Euro</u>
313005	0,05mg	
313025	0,25mg	
313100	1mg	

Version 1.1 07.03

Bioron GmbH

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916
E-Mail: info@bioron.net NET: www.bioron.net