

**Just fine Molecular Biology**

**DEFINITION** Cloning vector pUC19 complete sequence.


**ACCESSION** M77789

**VERSION** M77789.2

```

1 ggcccaata cgcaaaccgc ctctccccgc gcgttgccg attcattaat gcagctggca
61 cgacaggttt cccgactgga aagcgggag tgagcgaac gcaattaatg tgagtagct
121 cactcattag gcaccccagg cttacactt tatgcttccg gctcgtatgt tgttggaat
181 tgtgagcggg taacaatttc acacagggaa cagctatgac catgattacg ccaagcttgc
241 atgctgcag gtcgactcta gaggatcccc ggtaccgag ctggaattca ctggccgtcg
301 tttacaacg tcgtgactgg gaaaaccctg gcgttaccga acttaatcgc ctgcagcac
361 atccccctt cgccagctgg cgtaatagcg aagaggcccc caccgatcgc ccttccaac
421 agttgcgag cctgaatggc gaatggcgc tgatgcggta ttttctctt acgcatctgt
481 gcggtatttc acaccgata tggtgactc tcagtacaat ctgctctgat gccgatagt
541 taagccagcc cgcacaccg ccaacaccg ctgacgcgcc ctgacgggt tigtctctcc
601 cggcatccg ttacagacaa gctgtgaccg tctccgggag ctgcatgtgt cagaggttt
661 caccgtcatc accgaaacgc gcgagacgaa agggcctcgt gatacgccta ttttatagg
721 ttaatgcat gataataatg gtttctaga cgtcaggtag cactttcgg ggaatgtgc
781 gcggaacccc tattgttta ttttctaaa tacattcaaa tatgtatccg ctcagagac
841 aataaccctg ataatgctt caataatatt gaaaaaggaa gagtatgagt attcaacatt
901 tccgtctgc ccttattccc tttttgcgg cattttgctt tctgtttt gctcaccag
961 aaacgctgt gaaagtaaaa gatgctgaag atcagttggg tgcacgagtg ggtacatcg
1021 aactggatct caacagcggg aagatccttg agagtttccg ccccgaagaa cgttttcaa
1081 tgatgagcac ttttaagtt ctgctatgtg gcgcggtatt atcccgtatt gacgccgggc
1141 aagagaact cggtcgccc atacactatt ctcagaatga ctggttgag tactaccag
1201 tcacagaaaa gcactctacg gatggcatga cagtaagaga attatgcagt gctgccataa
1261 ccatgagtga taacactcgc gccaaactac ttctgacaac gatcggagga ccgaaggagc
1321 taaccgcttt ttgcacaac atgggggatc atgtaactcg cttgatcgt tgggaaccgg
1381 agctaatga agccatacca aacgacgagc gtacaccac gatgcctgta gcaatggcaa
1441 caacgttgcg caaactatta actggcgaac tactactct agcttcccg caacaattaa
1501 tagactggat ggagggcggat aaagttgcag gaccactct gcgctggcc ctccggctg
1561 gctggtttat tctgataaaa tctggagccg tgagcgtgg gctcgcggt atcattgcag
1621 cactggggcc agatgtaag ccttcccgtg tcgtagttat ctacacgacg gggagtcagg
1681 caactatgga tgaacgaaat agacagatcg ctgagatagg tgctcactg attaagcatt
1741 ggtaactgic agaccaagtt tactcatata tactttagat tgatttaaaa ctctatttt
1801 aatttaaag gatctagggt aagatccttt tgataatct catgacaaa atcccitac
1861 gtgagtttc gtccactga cgtcagacc ccgtagaaaa gatcaaagga tcttctgag
1921 atccttttt tctgcgcta atctgtgct tgcaaaaaa aaaaccaccg ctaccagcgg
1981 tggttgttt gccgatcaa gagctacca ctttttcc gaagtaact ggctcagca
2041 gagcgcagat accaaatact gtcttctag ttagccgta gtaggccac cactcaaga
2101 actctgtagc accgcctaca tacctcgtc tgctaactct gttaccagtg gctgctgcca
2161 gtggcgataa gctgtgctt accgggttg actcaagacg atagttaccg gataaggcgc
2221 agcggtcggg ctgaacgggg ggttcgtgca cacagcccag ctggagcga acgacctaca
2281 ccgaactgag atacctacag cgtgagctat gagaagcgc cagcttccc gaaggagaa

```



```

2341 agcggcagc gtagccgta agcggcagg tcggaacagg
agagcgcagc agggagctc
2401 caggggaaa cgctggtat cttatagtc ctgtcgggt
tcgccacct tgactgagc
2461 gtcgatttt gtgatgctg tcaggggggc ggagcctatg
gaaaacgcc agcaacgagg
2521 ccttttacg gttcctggcc tttgctggc ctttgctca catgttctt
cctgcgttat
2581 cccctgattc tgtggataac cgtattaccg cctttgagt
agctgatacc gctcggcga
2641 gccgaacgac cgagcgcagc gagtcagtga gcgaggaaga
ggaaga

```

**Bioron GmbH**

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916  
 E-Mail: [info@bioron.net](mailto:info@bioron.net) NET: [www.bioron.net](http://www.bioron.net)