



Just fine Molecular Biology

pBR322/HaeIII Marker DNA,

Description pBR322 digest of pBR 322 plasmid yields 22 fragments.

Fragments sizes (bp) 587, 540, 502, 458, 434, 267, 234, 213, 192, 184, 124, 123, 104, 89, 64, 57, 51, 21, 18, 11, 8


Storage buffer 10mM Tris-HCl (pH 7.8), 10mM NaCl, 1mM EDTA.

Concentration: 0,2 mg/ml

Storage at -20°C

Loading the recommended amount of size marker to load on an agarose gel is 1 µg/lane.

1,7% agarose



<u>Cat#</u>	<u>Pack size</u>	<u>Price, Euro</u>
302005	50 µg	
302025	250µg	

Version 1.1 12.04

Bioron GmbH

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916

E-Mail: info@bioron.net

NET: www.bioron.net



Just fine Molecular Biology

Bioron GmbH

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916

E-Mail: info@bioron.net

NET: www.bioron.net