



The ENZYME Company

Random Primers 30 µg

Description:

The primer is a mixture of single-stranded random hexanucleotides with 5'– and 3'–hydroxyl ends. Random Hexamer Primers can be used for first-strand cDNA synthesis and cloning. The primer is supplied in lyophilized form.

| Catalog # | Pack size | Price, Euro | SHOP |
|-----------|-----------|-------------|-----------------------------|
| 402001 | 30 µg | 35.00 | Add to cart |

Version 28.08.06 AS

Bioron GmbH

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916
E-Mail: info@bioron.net NET: www.bioron.net