



The ENZYME Company

« Complete II KCl » - Reaction buffer for :
Taq DNA polymerase

Reaction buffer (x10) complete II KCl:

500 mM KCl, 100 mM TrisHCl pH8,8, 0,1% Tween-20, 15 mM MgCl₂

Catalog #	Pack size
102006	1 ml
102007	5 ml

« Incomplete » - Reaction buffer for :

Bioron

Contact Germany Phone: +49-(0)-621- 5720 915

Contact Singapore Phone +65 6896 6942

E-Mail: info@bioron.net

WEB: www.bioron.net



The ENZYME Company

DFS-Taq DNA polymerase

Cat.-No. : 102002

5 x 1 ml

Reaction buffer (x10) incomplete:

160 mM (NH₄)₂SO₄, 670mM TrisHCl pH8,8, 0,1% Tween-20

Catalog #	Pack size
102001	1 ml
102002	5 ml

Bioron

Contact Germany Phone: +49-(0)-621- 5720 915

Contact Singapore Phone +65 6896 6942

E-Mail: info@bioron.net

WEB: www.bioron.net