



***The ENZYME Company***

---

« Incomplete » - Reaction buffer for :  
Taq DNA polymerase

*Reaction buffer (x10) incomplete:*

160 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 670mM TrisHCl pH8,8, 0,1% Tween-20

<b>Catalog #</b>	<b>Pack size</b>
102001	1 ml
102002	5 ml

---

« Incomplete » - Reaction buffer for :

---

**Bioron**

Contact Germany Phone: +49-(0)-621- 5720 915

Contact Singapore Phone +65 6896 6942

**E-Mail: [info@bioron.net](mailto:info@bioron.net)**

**WEB: [www.bioron.net](http://www.bioron.net)**



***The ENZYME Company***

**DFS-Taq DNA polymerase**

**Cat.-No. : 102002**

**5 x 1 ml**

***Reaction buffer (x10) incomplete:***

160 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 670mM TrisHCl pH8,8, 0,1% Tween-20

<b>Catalog #</b>	<b>Pack size</b>
102001	1 ml
102002	5 ml

---

**Bioron**

Contact Germany Phone: +49-(0)-621- 5720 915

Contact Singapore Phone +65 6896 6942

**E-Mail: [info@bioron.net](mailto:info@bioron.net)**

**WEB: [www.bioron.net](http://www.bioron.net)**