



The ENZYME Company

Enzyme	Prototype	Recognition Sequence	Cat No	Package, u.a.
--------	-----------	----------------------	--------	---------------

Bme18 I	Ava II	G <sup>^</sup> GWCC CCWG <sup>^</sup> G	266010	1000 units
			260050	5000 units

Lot-number:	Assayed:	Quantity:

Origin	Bacillus megaterium 18
Concentration	5000-10000, u.a./ml
Storage conditions	10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 ug/ml BSA; 50% glycerol; Store at -20°C.
Ligation	After 10-fold overdigestion with enzyme more than 90% of the DNA fragments can be ligated and recut.
Non-specific activity	No nonspecific activity was detected after incubation of 1 ug of DNA with 20 u.a. of enzyme for 16 hours at 37°C.
Optimum temperature	37 °C
Inactivation 20 minutes under 65 °C	Yes
Optimum SE-buffer	O (50 mM Tris-HCl (pH 7.6 at 25°C); 10 mM MgCl <sub>2</sub> ; 100 mM NaCl; 1 mM DTT.)

Enzyme activity in % of maximum :

B	G	O	W	Y
10 - 25	25 - 50	100	75 - 100	10 - 25

References: Degtyarev, S.Kh., Rechkunova, N.I., Grinev, A.A., Dedkov, V.S. Izv. Sib. Otd. Akad. Nauk SSSR 15: 25-26 (1989).

Unit-definition	One unit of the enzyme is the amount required to hydrolyze 1 µg of DNA in 1 hour in a total reaction volume of 50 µl. Concentrated enzymes are diluted to approximately 1000 units/ml with the buffer (10mM Tris-HCL (ph7.6); 50 mM KCL; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol) before determining their activity.
-----------------	---

**Bioron GmbH**

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916  
E-Mail: [info@bioron.net](mailto:info@bioron.net) NET: [www.bioron.net](http://www.bioron.net)