



The ENZYME Company

Enzyme	Prototype	Recognition Sequence	Cat No	Package, u.a.
--------	-----------	----------------------	--------	---------------

Bsp19 I	Nco I	C [^] CATGG GGTAC [^] C	215005	500 units
			215025	2500 units

Lot-number:	Assayed:	Quantity:
S18	02.2008	2500 units

Origin	Bacillus species 19
Concentration	20000 units/ml
Storage conditions	10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 ug/ml BSA; 50% glycerol. Store at -20°C.
Ligation	After 10-fold overdigestion with enzyme more than 90% of the DNA fragments can be ligated and recut.
Non-specific activity	No nonspecific activity was detected after incubation of 1 ug of DNA with 20 u.a. of enzyme for 16 hours at 37°C.
Optimum temperature	37 °C
Inactivation 20 minutes under 65 °C	Yes
Optimum SE-buffer	2W (20 mM Tris-HCl (pH 8,5 at 25°C); 10 mM MgCl ₂ ; 200 mM NaCl; 1 mM DTT + 100 ug/ml BSA) + BSA

Enzyme activity in % of maximum :

B	G	O	W	Y
0 - 10	10 - 25	50 - 75	75 - 100	10 - 25

Note: To obtain 100% activity, BSA should be added to the 1x reaction mix to a final concentration of 100 ug/ml.

References: Repin, V.E., Rechkunova, N.I., Degtyarev, S.Kh. Unpublished observations. Repin, V.E., Lebedev, L.R., Andreeva, I.S., Puchkova, L.I., Zernov, Y.P., Serov, G.D., Tereshchenko, T.A., Aphinogenova, G.N., Pustoshilova, N.M. Biotekhnologiya 0: 18-27 (1998).

Unit-definition	One unit of the enzyme is the amount required to hydrolyze 1 µg of DNA in 1 hour in a total reaction volume of 50 µl. Concentrated enzymes are diluted to approximately 1000 units/ml with the buffer (10mM Tris-HCL (ph7.6); 50 mM KCL; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol) before determining their activity.
------------------------	---

Bioron GmbH

Contact: Phone: +49-(0)-621- 5720 915 Fax:+49-(0)-621-5720 916
E-Mail: info@bioron.net NET: www.bioron.net